

NSK

NEW PRODUCT



Varios Combi Pro2

The NSK Varios Combi Pro2 (VCP2) is a powerful dual-function system designed for seamless oral hygiene and periodontal care. Combining supragingival and subgingival powder therapy with precision ultrasonic scaling, the NSK VCP2 is engineered for modern, minimally invasive workflows that help optimise clinical efficiency without compromising control or patient comfort.

VARIOS COMBI PRO2



INCLUDES 3 TIPS
(P20, G8, G16)



Basic Set

Y1500739

RRP: \$10,910

Includes:

- Varios Combi Pro2 Unit
- Prophy Air Polishing Kit (Supra)
- Wireless Foot Control
- CLM Biofilm Prevention Bottle Set

Accessories:

- Optic Titanium Scaler Handpiece
- Ultrasonic Scaler Tips (P20, G8, G16)
- Prophy Powder Handpiece
- 1x FLASH pearl powder bottle (Supra)



Bundle Package

VCP2BUNDLEPKG

RRP: \$17,599

Includes:

- Varios Combi Pro2 Basic Set
- Perio Air Polishing Kit (Sub) + 1x Jet Nozzle

Extra Accessories:

- 1x Optic Titanium Scaler Handpiece
- 2x Prophy Powder Handpieces
- 1x Perio Powder Handpiece
- 4x SOFT Pearl Powder
- 2x Perio Mate Powder
- 1x CLM Module Set
- 1x Water Bottle Kit
- 1x VCP Cart



DISCOVER AIR-POLISHING POWDERS



POWDERS

10% OFF RRP

FLASH pearl

Y900698

Calcium carbonate based powder (54µm) for stain removal.

4 x 300g bottles



RRP: \$343

NOW: \$309

SOFT pearl

Y1500721

Glycine based powder (65µm), ideal for moderate staining.

4 x 250g bottles



RRP: \$343

NOW: \$309

Perio Mate Powder

Y900938

Glycine based powder (25µm), for sub and supragingival use.

2 x 160g bottles



RRP: \$329

NOW: \$296

DISCOVER MORE TIPS



VARIOS TIPS

10% OFF RRP

SCALING TIPS

G1 Scaling Tip

Z217101

RRP: \$146

NOW: \$131



G2 Scaling Tip

Z217102

RRP: \$146

NOW: \$131



G3 Scaling Tip

Z217103

RRP: \$146

NOW: \$131



G4 Scaling Tip

Z217104

RRP: \$146

NOW: \$131



G5 Scaling Tip

Z217105

RRP: \$146

NOW: \$131



G6 Scaling Tip

Z217106

RRP: \$169

NOW: \$152



G8 Scaling Tip

Z217108

RRP: \$146

NOW: \$131




G16 Scaling Tip

Z217120

RRP: \$169

NOW: \$152



PERIO TIPS

P1 Perio Tip

Z217404

RRP: \$150

NOW: \$135



P2D Perio Tip

Z217402

RRP: \$157

NOW: \$141



P20 Perio Tip

Z217420

RRP: \$225

NOW: \$202



P10 Perio Tip

Z217410

RRP: \$225

NOW: \$202



P11R Perio Tip

Z217408

RRP: \$225

NOW: \$202



P11L Perio Tip

Z217409

RRP: \$225

NOW: \$202



P25L Perio Tip

Z217425

RRP: \$225

NOW: \$202



P25R Perio Tip

Z217424

RRP: \$225

NOW: \$202

